

For Research Use Only

GOLGA7 Polyclonal antibody

Catalog Number: 14723-1-AP

Featured Product



Basic Information

Catalog Number:

14723-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6460

GenBank Accession Number:

BC012032

GeneID (NCBI):

51125

UNIPROT ID:

Q7Z5G4

Full Name:

golgi autoantigen, golgin subfamily a, 7

Calculated MW:

16 kDa

Observed MW:

15-20 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

IHC 1:50-1:500

Applications

Tested Applications:

IHC, WB, ELISA

Species Specificity:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : A549 cells, HEK-293 cells, HL-60 cells, HepG2 cells

IHC : mouse stomach tissue,

Background Information

GOLGA7 (Golgin subfamily A member 7) is an acylated Golgi protein that interacts with GCP170 protein (PMID: 14522980). The ZDHHC9-GOLGA7 complex is a palmitoyltransferase specific for HRAS and NRAS. The ZDHHC5-GOLGA7 protein acyltransferase complex promotes nonapoptotic cell death (PMID: 31631010).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 14723-1-AP (GOLGA7 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.



WB result of GOLGA7 antibody (14723-1-AP; 1:400; incubated at room temperature for 1.5 hours) with sh-Control and sh-GOLGA7 transfected HEK-293 cells.



Immunohistochemical analysis of paraffin-embedded mouse stomach tissue slide using 14723-1-AP (GOLGA7 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse stomach tissue slide using 14723-1-AP (GOLGA7 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).