

For Research Use Only

PAFAH1B2 Polyclonal antibody

Catalog Number: 14729-1-AP

Featured Product

2 Publications



Basic Information

Catalog Number:

14729-1-AP

Size:

180 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6489

GenBank Accession Number:

BC000398

GeneID (NCBI):

5049

UNIPROT ID:

P68402

Full Name:

platelet-activating factor
acetylhydrolase, isoform Ib, beta
subunit 30kDa

Calculated MW:

30 kDa

Observed MW:

30 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

IP 0.5-4.0 µg for 1.0-3.0 mg of total
protein lysate

Applications

Tested Applications:

IP, WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse

Cited Species:

human, mouse

Positive Controls:

WB: human testis tissue, mouse small intestine tissue

IP: mouse testis tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
YingYing He	34930255	Cancer Cell Int	WB,IHC
Yasuhiro Nakamura	36724953	Biol Pharm Bull	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 14729-1-AP (PAFAH1B2 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.



IP result of anti-PAFAH1B2 (IP:14729-1-AP, 4ug; Detection:14729-1-AP 1:600) with mouse testis tissue lysate 4000ug.