

For Research Use Only

TWF2 Polyclonal antibody

Catalog Number: 14769-1-AP

1 Publications



Basic Information

Catalog Number:

14769-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6553

GenBank Accession Number:

BC000327

GeneID (NCBI):

11344

UNIPROT ID:

Q6IBS0

Full Name:

twinfilin, actin-binding protein,
homolog 2 (Drosophila)

Calculated MW:

40 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Zhang Kaishang	30391462	Life Sci	WB, IHC

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

