

For Research Use Only

NOL4 Polyclonal antibody

Catalog Number: 14802-1-AP **1 Publications**



Basic Information

Catalog Number:

14802-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6518

GenBank Accession Number:

BC000313

GeneID (NCBI):

8715

UNIPROT ID:

O94818

Full Name:

nucleolar protein 4

Calculated MW:

71 kDa

Observed MW:

64-70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

rat

Positive Controls:

WB : mouse brain tissue, mouse testis tissue, rat brain tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Masaki Hiramoto	29207083	Int J Mol Med	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 14802-1-AP (NOL4 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.