

For Research Use Only

UBOX5 Polyclonal antibody

Catalog Number: 14815-1-AP



Basic Information

Catalog Number:

14815-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6565

GenBank Accession Number:

BC000515

GeneID (NCBI):

22888

UNIPROT ID:

O94941

Full Name:

U-box domain containing 5

Calculated MW:

59 kDa

Observed MW:

59-69 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse kidney tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 14815-1-AP (UBO \times 5 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.