

For Research Use Only

# TAB1 Polyclonal antibody

Catalog Number: 14819-1-AP

Featured Product

5 Publications



## Basic Information

Catalog Number:

14819-1-AP

Size:

380 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6590

GenBank Accession Number:

BC050554

GeneID (NCBI):

10454

UNIPROT ID:

Q15750

Full Name:

mitogen-activated protein kinase  
kinase kinase 7 interacting protein 1

Calculated MW:

55 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, ColP

Species Specificity:

human, mouse, rat

Cited Species:

human, sheep

## Background Information

### Notable Publications

| Author         | Pubmed ID | Journal              | Application |
|----------------|-----------|----------------------|-------------|
| Zhongjia Jiang | 28553017  | Mediators Inflamm    | WB          |
| Lin Yuan       | 35352799  | Cardiovasc Res       | WB, ColP    |
| Kun-Lin Chen   | 31285783  | Oxid Med Cell Longev | WB          |

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

