

For Research Use Only

HMHA1 Polyclonal antibody

Catalog Number: 14832-1-AP

1 Publications



Basic Information

Catalog Number:

14832-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6636

GenBank Accession Number:

BC065223

GeneID (NCBI):

23526

UNIPROT ID:

Q92619

Full Name:

histocompatibility (minor) HA-1

Calculated MW:

1136 aa, 125 kDa

Observed MW:

140-150 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

IF

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Hongyu Ma	37199734	Aging Dis	IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

