

For Research Use Only

HMHA1 Polyclonal antibody, PBS Only

Catalog Number:14832-1-PBS



Basic Information

Catalog Number:

14832-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6636

GenBank Accession Number:

BC065223

GeneID (NCBI):

23526

UNIPROT ID:

Q92619

Full Name:

histocompatibility (minor) HA-1

Calculated MW:

1136 aa, 125 kDa

Observed MW:

140-150 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data