

For Research Use Only

ZIC4-Specific Polyclonal antibody

Catalog Number: 14835-1-AP



Basic Information

Catalog Number:

14835-1-AP

Size:

90 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_032153

GeneID (NCBI):

84107

UNIPROT ID:

Q8N9L1

Full Name:

Zic family member 4

Calculated MW:

37 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

The ZIC gene family includes 5 genes that are highly conserved across evolution. These genes encode zinc finger proteins that are expressed during development and maturation of the central nervous system and have critical roles in the development of the cerebellum. ZIC4, like other members, may have roles in the nervous system development.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

