

For Research Use Only

GTF2H3 Polyclonal antibody

Catalog Number: 14857-1-AP



Basic Information

Catalog Number:

14857-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6652

GenBank Accession Number:

BC065250

GeneID (NCBI):

2967

UNIPROT ID:

Q13889

Full Name:

general transcription factor IIH,
polypeptide 3, 34kDa

Calculated MW:

34 kDa

Observed MW:

30-34 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HL-60 cells, HT-29 cells, HeLa cells, HepG2 cells,
MCF-7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 14857-1-AP (GTF2H3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.