

For Research Use Only

URAT1 Polyclonal antibody, PBS Only

Catalog Number: 14937-1-PBS

Featured Product



Basic Information

Catalog Number:

14937-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6770

GenBank Accession Number:

BC053348

GeneID (NCBI):

116085

UNIPROT ID:

Q96S37

Full Name:

solute carrier family 22 (organic anion/urate transporter), member 12

Calculated MW:

60 kDa

Observed MW:

65-70 kDa, 40 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat, pig

Background Information

Human urate anion exchanger (URAT1), also known as SLC22A12, is a urate transporter located in the proximal tubules of human kidney and responsible for urate reabsorption. Defects in gene encoding URAT1 cause hypouricemia renal type 1 (RHUC1). This antibody recognizes the native and glycosylated URAT1 between 65 and 70 kDa as well as the truncated 40kDa protein.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 14937-1-AP (URAT1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 14937-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using 14937-1-AP (URAT1 antibody) at dilution of 1:400 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 14937-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using 14937-1-AP (URAT1 antibody) at dilution of 1:400 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 14937-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 14937-1-AP (URAT1 antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 14937-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 14937-1-AP (URAT1 antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 14937-1-PBS in a different storage buffer formulation.