

For Research Use Only

# SMAGP Polyclonal antibody

Catalog Number: 14998-1-AP



## Basic Information

Catalog Number:

14998-1-AP

Size:

180 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6988

GenBank Accession Number:

BC003379

GeneID (NCBI):

57228

UNIPROT ID:

Q0VAQ4

Full Name:

small trans-membrane and glycosylated protein

Calculated MW:

6 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

Small cell transmembrane and glycosylated protein (SMAGP) is an evolutionary conserved protein expressed on the lateral epithelial cell membrane. Altered expression of SMAGP has been reported in cancer, suggesting that SMAGP may have a role in tumor progression (PMID: 15021913; 15986429). This polyclonal antibody is raised against LOC57228 protein (GenBank: AAC72956.1; Uniprot Database Entry: O95332).

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

