

For Research Use Only

# UNKL Polyclonal antibody, PBS Only

Catalog Number:15000-1-PBS



## Basic Information

Catalog Number:

15000-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6990

GenBank Accession Number:

BC000150

GeneID (NCBI):

64718

UNIPROT ID:

Q9H9P5

Full Name:

unkempt homolog (Drosophila)-like

Calculated MW:

25 kDa

Observed MW:

32 kDa, 45 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

Putative E3 ubiquitin-protein ligase UNKL was expressed at significantly lower levels in endometrial tumor tissues as compared to the endometrium. Importantly, in human endometrial cancer, primary tumor expression of UNKL was correlated with overall survival in black patients with high mutational burden.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 15000-1-AP (UNKL antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 15000-1-PBS in a different storage buffer formulation.