

For Research Use Only

USP10 Polyclonal antibody

Catalog Number: 15025-1-AP



Basic Information

Catalog Number: 15025-1-AP	GenBank Accession Number: BC000263	Purification Method: Antigen affinity purification
Size: 490 µg/ml	GeneID (NCBI): 9100	
Source: Rabbit	UNIPROT ID: Q14694	
Isotype: IgG	Full Name: ubiquitin specific peptidase 10	
Immunogen Catalog Number: AG7106	Calculated MW: 87 kDa, 92 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

USP10, also named as USP32 and NY-REN-60, belongs to the peptidase C19 family. The deubiquitinating protease USP10 has been identified as a new regulator of p53 in the DNA damage response and tumor development.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data