

For Research Use Only

dNGLUC Polyclonal antibody

Catalog Number: 15028-1-AP



Basic Information

Catalog Number:

15028-1-AP

Size:

130 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7427

GenBank Accession Number:

AY015993

GeneID (NCBI):

Full Name:

Calculated MW:

20 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

Background Information

Luciferase is a generic term for the class of oxidative enzymes used in bioluminescence and is distinct from a photoprotein. The immunogen of this antibody the luciferase (AY015993) from the copepod marine organism *Gaussia princeps*. This antibody is used to test the luciferase.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

