

For Research Use Only

CUTC Polyclonal antibody, PBS Only

Catalog Number:15030-1-PBS



Basic Information

Catalog Number:

15030-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2126

GenBank Accession Number:

BC021105

GeneID (NCBI):

51076

UNIPROT ID:

Q9NTM9

Full Name:

cutC copper transporter homolog (E. coli)

Calculated MW:

29 kDa

Observed MW:

29 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




HeLa cells were subjected to SDS PAGE followed by western blot with 15030-1-AP (CUTC antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 15030-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human endometrial cancer using 15030-1-AP (CUTC antibody) at dilution of 1:100 (under 10x lens). This data was developed using the same antibody clone with 15030-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human endometrial cancer using 15030-1-AP (CUTC antibody) at dilution of 1:100 (under 40x lens). This data was developed using the same antibody clone with 15030-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CUTC antibody (15030-1-AP) at dilution of 1:400 and Coralite@488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-phalloidin (red). This data was developed using the same antibody clone with 15030-1-PBS in a different storage buffer formulation.