

For Research Use Only

EMP1 Polyclonal antibody, PBS Only

Catalog Number:15038-1-PBS



Basic Information

Catalog Number:

15038-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6771

GenBank Accession Number:

BC047300

GeneID (NCBI):

2012

UNIPROT ID:

P54849

Full Name:

epithelial membrane protein 1

Calculated MW:

18 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data