

For Research Use Only

DYNLRB1 Polyclonal antibody

Catalog Number: 15049-1-AP **3 Publications**



Basic Information

Catalog Number: 15049-1-AP	GenBank Accession Number: BC002481	Purification Method: Antigen affinity purification
Size: 200 µg/ml	GeneID (NCBI): 83658	
Source: Rabbit	UNIPROT ID: Q9NP97	
Isotype: IgG	Full Name: dynein, light chain, roadblock-type 1	
Immunogen Catalog Number: AG7032	Calculated MW: 11 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
IF

Species Specificity:
human, mouse, rat

Cited Species:
human, rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Syeda Ridita Sharif	27646688	Mol Cells	IF
Md Ariful Islam	26467288	Mol Cells	IF
Md Ariful Islam	26272270	Exp Mol Med	IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data