

For Research Use Only

CINP Polyclonal antibody

Catalog Number: 15050-1-AP



Basic Information

Catalog Number:

15050-1-AP

Size:

670 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7033

GenBank Accession Number:

BC000600

GeneID (NCBI):

51550

UNIPROT ID:

Q9BW66

Full Name:

cyclin-dependent kinase 2-interacting protein

Calculated MW:

24 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

CINP is a cell-cycle checkpoint protein. It interacts with ATR interacting protein and regulates ATR-dependent signaling, resistance to replication stress, and G2 checkpoint integrity.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

