

For Research Use Only

DHRS7 Polyclonal antibody, PBS Only

Catalog Number:15061-1-PBS



Basic Information

Catalog Number:

15061-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7064

GenBank Accession Number:

BC000637

GeneID (NCBI):

51635

UNIPROT ID:

Q9Y394

Full Name:

dehydrogenase/reductase (SDR family) member 7

Calculated MW:

38 kDa

Observed MW:

32-38 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

DHRS7 (Dehydrogenase/Reductase (SDR family) member 7), also known as SDR34C1, retSDR4, is a member of a large family of short-chain dehydrogenases/reductases (SDR). The DHRS7 gene codes for two subtypes and is found on chromosome 14. Isotype 1 (38 kDa) has 339 amino acids, while isotype 2 (32 kDa) contains 289 amino acids (PMID: 26466768). DHRS7 has been demonstrated to be an integral protein facing the lumen of the endoplasmic reticulum with lack of posttranscriptional glycosylation modification (PMID: 24246760).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



LNCaP cells were subjected to SDS PAGE followed by western blot with 15061-1-AP (DHR57 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 15061-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded rat skeletal muscle tissue slide using 15061-1-AP (DHR57 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 15061-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse skeletal muscle tissue slide using 15061-1-AP (DHR57 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 15061-1-PBS in a different storage buffer formulation.