

For Research Use Only

APC11 Polyclonal antibody

Catalog Number: 15065-1-AP



Basic Information

Catalog Number:

15065-1-AP

Size:

470 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7069

GenBank Accession Number:

BC000607

GeneID (NCBI):

51529

UNIPROT ID:

Q9NYG5

Full Name:

anaphase promoting complex subunit

11

Calculated MW:

9 kDa, 21 kDa

Observed MW:

21 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse pancreas tissue, A549 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with 15065-1-AP (APC11 antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.