

For Research Use Only

LPCAT2 Polyclonal antibody

Catalog Number: 15082-1-AP **1 Publications**



Basic Information

Catalog Number:

15082-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7156

GenBank Accession Number:

BC002472

GeneID (NCBI):

54947

UNIPROT ID:

Q7L5N7

Full Name:

lysophosphatidylcholine
acyltransferase 2

Calculated MW:

60 kDa

Observed MW:

53 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:5000-1:50000

IF 1:200-1:800

Applications

Tested Applications:

IF/ICC, WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : HEK-293 cells, LPS treated HEK-293, HeLa cells,
Jurkat cells, RAW 264.7 cells

IF : A431 cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Melanie Merkel	37627584	Antioxidants (Basel)	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 15082-1-AP (LPCAT2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Ethanol) fixed A431 cells using LPCAT2 antibody (15082-1-AP) at dilution of 1:400 and CoraLite® 488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L), CL594-Phalloidin (red).