

For Research Use Only

ALG12 Polyclonal antibody

Catalog Number: 15155-1-AP



Basic Information

Catalog Number:

15155-1-AP

Size:

430 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7201

GenBank Accession Number:

BC001729

GeneID (NCBI):

79087

UNIPROT ID:

Q9BV10

Full Name:

asparagine-linked glycosylation 12,
alpha-1,6-mannosyltransferase
homolog (S. cerevisiae)

Calculated MW:

55 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

