

For Research Use Only

# GALK2 Polyclonal antibody

Catalog Number: 15164-1-AP

Featured Product



## Basic Information

Catalog Number:

15164-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7278

GenBank Accession Number:

BC005141

GeneID (NCBI):

2585

UNIPROT ID:

Q01415

Full Name:

galactokinase 2

Calculated MW:

50 kDa

Observed MW:

45-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HL-60 cells, HEK-293 cells, HeLa cells, HepG2 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HL-60 cells were subjected to SDS PAGE followed by western blot with 15164-1-AP (GALK2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



WB result of GALK2 antibody (15164-1-AP; 1:2000; incubated at room temperature for 1.5 hours) with sh-Control and sh-GALK2 transfected HepG2 cells.