

For Research Use Only

GEMIN8 Polyclonal antibody, PBS Only

Catalog Number:15251-1-PBS



Basic Information

Catalog Number:

15251-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6774

GenBank Accession Number:

BC003607

GeneID (NCBI):

54960

UNIPROT ID:

Q9NWZ8

Full Name:

gem (nuclear organelle) associated protein 8

Calculated MW:

29 kDa

Observed MW:

32 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IF/ICC, IP, Indirect ELISA

Species Specificity:

human

Background Information

GEMIN8, also named as FAM51A1 and Gemin-8, is associated with the SMN complex and is localized in the cytoplasm and in the nucleus, where it is highly concentrated in Gems. Gemin6, Gemin7, and Unrip form a stable cytoplasmic complex whose association with SMN requires Gemin8. Gemin8 binds directly to SMN and mediates its interaction with the Gemin6/Gemin7 heterodimer. This antibody recognize GEMIN8(28-32kd) only.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human liver cancer using 15251-1-AP (GEMIN8 antibody) at dilution of 1:50 (under 10x lens). This data was developed using the same antibody clone with 15251-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human liver cancer using 15251-1-AP (GEMIN8 antibody) at dilution of 1:50 (under 40x lens). This data was developed using the same antibody clone with 15251-1-PBS in a different storage buffer formulation.



IP result of anti-GEMIN8 (IP:15251-1-AP, 3ug; Detection:15251-1-AP 1:400) with HepG2 cells lysate 1720ug. This data was developed using the same antibody clone with 15251-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using 15251-1-AP (GEMIN8 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). This data was developed using the same antibody clone with 15251-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 15251-1-AP (GEMIN8 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 15251-1-PBS in a different storage buffer formulation.