

TRIM62 Polyclonal antibody

Catalog Number: 15256-1-AP

1 Publications

Basic Information

Catalog Number:

15256-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7143

GenBank Accession Number:

BC001222

GeneID (NCBI):

55223

UNIPROT ID:

Q9BVG3

Full Name:

tripartite motif-containing 62

Calculated MW:

54 kDa

Observed MW:

54 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

IHC 1:50-1:500

IF/ICC 1:200-1:800

Applications

Tested Applications:

WB, IHC, IF/ICC, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : BGC-823 cells, SGC-7901 cells, mouse ovary tissue

IHC : mouse liver tissue, rat ovary tissue

IF/ICC : U2OS cells,

Background Information

TRIM62(Tripartite motif-containing protein 62) belongs to the TRIM/RBCC family. TRIM proteins are involved in a broad range of biological processes, including some that underlie genetic disorders, neurological disorders and cancers. TRIM62 is shown to promote the release of HIV from infected cells, and a TRIM62 protein that has a mutation in its RING domain was more active than the native protein in promoting HIV release(PMID:18836477).

Notable Publications

Author	Pubmed ID	Journal	Application
Uchil Pradeep D PD	23077300	J Virol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 4006900926


E: Proteintech-CN@ptglab.comW: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 15256-1-AP (TRIM62 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.




Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using 15256-1-AP (TRIM62 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).




Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using 15256-1-AP (TRIM62 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded rat ovary tissue slide using 15256-1-AP (TRIM62 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded rat ovary tissue slide using 15256-1-AP (TRIM62 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (-20°C Ethanol) fixed U2OS cells using TRIM62 antibody (15256-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-phalloidin (red).