

For Research Use Only

SMPD2 Polyclonal antibody

Catalog Number: 15290-1-AP

1 Publications



Basic Information

Catalog Number:

15290-1-AP

Size:

107 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7482

GenBank Accession Number:

BC000038

GeneID (NCBI):

6610

UNIPROT ID:

O60906

Full Name:

sphingomyelin phosphodiesterase 2,
neutral membrane (neutral
sphingomyelinase)

Calculated MW:

48 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

IF

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Wei Fan	34042046	Elife	IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

