

For Research Use Only

RRS1 Polyclonal antibody, PBS Only

Catalog Number: 15329-1-PBS

Featured Product



Basic Information

Catalog Number:

15329-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7030

GenBank Accession Number:

BC001811

GeneID (NCBI):

23212

UNIPROT ID:

Q15050

Full Name:

RRS1 ribosome biogenesis regulator homolog (*S. cerevisiae*)

Calculated MW:

41 kDa

Observed MW:

41 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



K-562 cells were subjected to SDS PAGE followed by western blot with 15329-1-AP (RRS1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 15329-1-PBS in a different storage buffer formulation.



IP result of anti-RRS1 (IP:15329-1-AP, 4ug; Detection:15329-1-AP 1:800) with K-562 cells lysate 3200ug. This data was developed using the same antibody clone with 15329-1-PBS in a different storage buffer formulation.



WB result of RRS1 antibody (15329-1-AP; 1:2000; incubated at room temperature for 1.5 hours) with sh-Control and sh-RRS1 transfected A549 cells. This data was developed using the same antibody clone with 15329-1-PBS in a different storage buffer formulation.