

For Research Use Only

# TAF1A Polyclonal antibody

Catalog Number: 15362-1-AP



## Basic Information

Catalog Number:

15362-1-AP

Size:

230 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8059

GenBank Accession Number:

BC013808

GeneID (NCBI):

9015

UNIPROT ID:

Q15573

Full Name:

TATA box binding protein (TBP)-  
associated factor, RNA polymerase I,  
A, 48kDa

Calculated MW:

53 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

## Background Information

TAF1A, also named as TAFI48 and SL1, is a component of the transcription factor SL1/TIF-IB complex. TAF1A has two isoform with MW 40 kDa and 48-52 kDa.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data