

For Research Use Only

LRRK1 Polyclonal antibody

Catalog Number: 15367-1-AP



Basic Information

Catalog Number:

15367-1-AP

Size:

250 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2599

GenBank Accession Number:

BC005408

GeneID (NCBI):

79705

UNIPROT ID:

Q385D2

Full Name:

leucine-rich repeat kinase 1

Calculated MW:

225 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data