

For Research Use Only

EEF1AKMT4 Polyclonal antibody

Catalog Number: 15418-1-AP

1 Publications



Basic Information

Catalog Number:

15418-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7706

GenBank Accession Number:

BC005835

GeneID (NCBI):

110599564

UNIPROT ID:

P0DPD7

Full Name:

EEF1A lysine methyltransferase 4

Calculated MW:

28 kDa

Observed MW:

28 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Positive Controls:

WB : mouse liver tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Mitchell Andrew J AJ	23894327	PLoS One	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 15418-1-AP (EEF1AKMT4 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.