

For Research Use Only

EEF1AKMT4 Polyclonal antibody

Catalog Number: 15418-1-AP **1 Publications**



Basic Information

Catalog Number: 15418-1-AP	GenBank Accession Number: BC005835	Purification Method: Antigen affinity purification
Size: 300 µg/ml	GeneID (NCBI): 110599564	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	UNIPROT ID: PODPD7	
Isotype: IgG	Full Name: EEF1A lysine methyltransferase 4	
Immunogen Catalog Number: AG7706	Calculated MW: 28 kDa	
	Observed MW: 28 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse liver tissue,
Cited Applications: WB	
Species Specificity: human, mouse, rat	
Cited Species: mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Mitchell Andrew J AJ	23894327	PLoS One	WB


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 15418-1-AP (EEF1AKMT4 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.