

For Research Use Only

HOXC8 Polyclonal antibody

Catalog Number: 15448-1-AP

Featured Product

9 Publications



Basic Information

Catalog Number:

15448-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7675

GenBank Accession Number:

BC053898

GeneID (NCBI):

3224

UNIPROT ID:

P31273

Full Name:

homeobox C8

Calculated MW:

27 aa, 7 kDa

Observed MW:

33-38 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF, IHC, CoIP

Species Specificity:

human

Cited Species:

human, sheep, mouse

Positive Controls:

WB : COLO 320 cells, MDA-MB-453s cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Masami Shiimori	34169237	iScience	IF
Joshua E Burda	35614216	Nature	IHC
Jie Zhang	31202850	Int J Biochem Cell Biol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



COLO 320 cells were subjected to SDS PAGE followed by western blot with 15448-1-AP (HOXC8 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.