

For Research Use Only

# HOXC8 Polyclonal antibody

Catalog Number: 15448-1-AP

Featured Product

8 Publications



## Basic Information

Catalog Number:

15448-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7675

GenBank Accession Number:

BC053898

GeneID (NCBI):

3224

UNIPROT ID:

P31273

Full Name:

homeobox C8

Calculated MW:

27 aa, 7 kDa

Observed MW:

33-38 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IF, IHC, CoIP

Species Specificity:

human

Cited Species:

human, sheep, mouse

Positive Controls:

WB : COLO 320 cells, MDA-MB-453s cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Masami Shiimori	34169237	iScience	IF
Joshua E Burda	35614216	Nature	IHC
Jie Zhang	31202850	Int J Biochem Cell Biol	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



COLO 320 cells were subjected to SDS PAGE followed by western blot with 15448-1-AP (HOXC8 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.