

For Research Use Only



# APOBEC3C-specific Polyclonal antibody

Catalog Number: 15495-1-AP

## Basic Information

Catalog Number:

15495-1-AP

Size:

310 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_014508

GeneID (NCBI):

27350

UNIPROT ID:

Q9NRW3

Full Name:

apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3C

Calculated MW:

23 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse

## Background Information

This antibody is specific to APOBEC3, not cross react with APOBEC1 and APOBEC2.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

