

For Research Use Only

EPHB6-Specific Polyclonal antibody

Catalog Number: 15498-1-AP



Basic Information

Catalog Number:

15498-1-AP

Size:

70 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_004445

GeneID (NCBI):

2051

UNIPROT ID:

O15197

Full Name:

EPH receptor B6

Calculated MW:

111 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

EPHB6 belongs to the protein kinase superfamily, Tyr protein kinase family and Ephrin receptor subfamily. EPHB6 modulates cell adhesion and migration by exerting both positive and negative effects upon stimulation with ephrin-B2. EPHB6 inhibits JNK activation, T cell receptor-induced IL-2 secretion and CD25 expression upon stimulation with ephrin-B2.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

