

For Research Use Only

HAGHL Polyclonal antibody, PBS Only

Catalog Number:15514-1-PBS



Basic Information

Catalog Number:

15514-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7729

GenBank Accession Number:

BC004353

GeneID (NCBI):

84264

UNIPROT ID:

Q6PII5

Full Name:

hydroxyacylglutathione hydrolase-like

Calculated MW:

32 kDa

Observed MW:

31 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

HAGHL, or HAGHL gene, is a gene located on chromosome 16p13.3, and it encodes a protein that has been identified as a novel metabolic oncogene involved in the progression of human colorectal cancer (CRC). Recent studies have demonstrated that HAGHL plays a critical role in regulating CRC by influencing various metabolic pathways, which may contribute to tumor growth and development. (PMID: 36417085)

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cells were subjected to SDS PAGE followed by western blot with 15514-1-AP (HAGHL antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 15514-1-PBS in a different storage buffer formulation.