

For Research Use Only

# SPTLC2 Polyclonal antibody

Catalog Number: 15515-1-AP

1 Publications



## Basic Information

Catalog Number:

15515-1-AP

Concentration:

120 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7774

GenBank Accession Number:

BC005123

GeneID (NCBI):

9517

UNIPROT ID:

O15270

Full Name:

serine palmitoyltransferase, long chain base subunit 2

Calculated MW:

63 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

## Background Information

## Notable Publications

Author	Pubmed ID	Journal	Application
Dongdong Huang	33411679	Aging (Albany NY)	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

