

For Research Use Only

# RPS25 Polyclonal antibody

Catalog Number: 15566-1-AP

Featured Product

1 Publications



## Basic Information

Catalog Number:

15566-1-AP

Size:

160 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7825

GenBank Accession Number:

BC004986

GeneID (NCBI):

6230

UNIPROT ID:

P62851

Full Name:

ribosomal protein S25

Calculated MW:

14 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Cited Applications:

WB, IP

Species Specificity:

human

Cited Species:

human

## Background Information

Ribosomal protein S25 (RPS25) is a novel MDM2 interacting protein. It may be involved in viral replication of Dicrostoviridae and hepatitis C viruses. RPS25 is overexpressed in human leukemia cells exhibiting adriamycin resistance. It may have a role in p53-mediated apoptosis and cell-cycle arrest. The expected MW of RPS25 is about 15-17 kDa. Catalog# 15566-1-AP is a rabbit polyclonal antibody raised against full-length human RPS25.

## Notable Publications

Author	Pubmed ID	Journal	Application
Zhang X X	22777350	Oncogene	WB, IP

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

