

For Research Use Only

PPP3CC Polyclonal antibody

Catalog Number: 15582-1-AP



Basic Information

Catalog Number:

15582-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8026

GenBank Accession Number:

BC004864

GeneID (NCBI):

5533

UNIPROT ID:

P48454

Full Name:

protein phosphatase 3 (formerly 2B), catalytic subunit, gamma isoform

Calculated MW:

57 kDa

Observed MW:

57 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB: Jurkat cells, K-562 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 15582-1-AP (PPP3CC antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.