

For Research Use Only

LSP1 Polyclonal antibody, PBS Only

Catalog Number:15593-1-PBS



Basic Information

Catalog Number:

15593-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7944

GenBank Accession Number:

BC001785

GeneID (NCBI):

4046

UNIPROT ID:

P33241

Full Name:

Lymphocyte-specific protein 1

Calculated MW:

37 kDa

Observed MW:

52 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, IP, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 15593-1-AP (LSP1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 15593-1-PBS in a different storage buffer formulation.



IP result of anti-LSP1 (IP:15593-1-AP, 4ug; Detection:15593-1-AP 1:300) with Raji cells lysate 4000ug. This data was developed using the same antibody clone with 15593-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 15593-1-AP (LSP1 antibody) at dilution of 1:800 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 15593-1-PBS in a different storage buffer formulation.