

For Research Use Only

CKS2 Polyclonal antibody

Catalog Number: 15616-1-AP **1 Publications**



Basic Information

Catalog Number: 15616-1-AP	GenBank Accession Number: BC006458	Purification Method: Antigen affinity purification
Size: 130 µg/ml	GeneID (NCBI): 1164	
Source: Rabbit	UNIPROT ID: P33552	
Isotype: IgG	Full Name: CDC28 protein kinase regulatory subunit 2	
Immunogen Catalog Number: AG8021	Calculated MW: 10 kDa	
	Observed MW: 10 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human, mouse, rat

Cited Species:
human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Long Jiang	28502657	Curr Biol	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data