

For Research Use Only

ICA1L Polyclonal antibody

Catalog Number: 15622-1-AP



Basic Information

Catalog Number:

15622-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8048

GenBank Accession Number:

BC000993

GeneID (NCBI):

130026

UNIPROT ID:

Q8NDH6

Full Name:

islet cell autoantigen 1,69kDa-like

Calculated MW:

54 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF/ICC 1:10-1:100

Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

IF/ICC : MCF-7 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of MCF-7 cells, using ICA1L antibody 15622-1-AP at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).