

For Research Use Only

# ICA1L Polyclonal antibody

Catalog Number: 15622-1-AP



## Basic Information

<b>Catalog Number:</b> 15622-1-AP	<b>GenBank Accession Number:</b> BC000993	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 500 µg/ml	<b>GeneID (NCBI):</b> 130026	<b>Recommended Dilutions:</b> IF 1:10-1:100
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q8NDH6	
<b>Isotype:</b> IgG	<b>Full Name:</b> islet cell autoantigen 1,69kDa-like	
<b>Immunogen Catalog Number:</b> AG8048	<b>Calculated MW:</b> 54 kDa	

## Applications

<b>Tested Applications:</b> IF/ICC,ELISA	<b>Positive Controls:</b> IF : MCF-7 cells,
<b>Species Specificity:</b> human, mouse, rat	

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of MCF-7 cells, using ICA1L antibody 15622-1-AP at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).