

For Research Use Only

IL-17RC Polyclonal antibody

Catalog Number: 15648-1-AP

2 Publications



Basic Information

Catalog Number:

15648-1-AP

Size:

147 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8190

GenBank Accession Number:

BC006411

GeneID (NCBI):

84818

UNIPROT ID:

Q8NAC3

Full Name:

interleukin 17 receptor C

Calculated MW:

791aa, 86 kDa; 538aa, 59 kDa

Observed MW:

45-50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human, mouse

Positive Controls:

WB : HEK-293 cells, HepG2 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xin Wu	35343786	mBio	WB
Panhong Wu	37086998	Sci Total Environ	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 15648-1-AP (IL17RC antibody) at dilution of 1:200 incubated at room temperature for 1.5 hours.