

For Research Use Only

PIGQ Polyclonal antibody

Catalog Number: 15658-1-AP



Basic Information

Catalog Number:

15658-1-AP

Size:

190 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8051

GenBank Accession Number:

BC006377

GeneID (NCBI):

9091

UNIPROT ID:

Q9BRB3

Full Name:

phosphatidylinositol glycan anchor
biosynthesis, class Q

Calculated MW:

84 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

Many eukaryotic proteins are anchored to the membrane through the hydrocarbon chains of a covalently bound glycosylphosphatidylinositol (GPI) membrane anchor. This gene encodes GPI1, also named as PIGQ, is involved in the first step of GPI biosynthesis. Multiple transcripts including calculated MW of 84kd, 65kd and 32kd were reported

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

