

For Research Use Only

# NOL6 Polyclonal antibody

Catalog Number: 15706-1-AP

1 Publications



## Basic Information

Catalog Number:

15706-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8378

GenBank Accession Number:

BC008298

GeneID (NCBI):

65083

UNIPROT ID:

Q9H6R4

Full Name:

nucleolar protein family 6 (RNA-associated)

Calculated MW:

1146 aa, 128 kDa

Observed MW:

130 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : SGC-7901 cells, HeLa cells, PC-3 cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Degang Dong	32249364	Cytotechnology	WB, IHC

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



SGC-7901 cells were subjected to SDS PAGE followed by western blot with 15706-1-AP (NOL6 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.