

For Research Use Only

ST6GALNAC6 Polyclonal antibody

Catalog Number: 15730-1-AP



Basic Information

Catalog Number:

15730-1-AP

Size:

210 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8475

GenBank Accession Number:

BC007802

GeneID (NCBI):

30815

UNIPROT ID:

Q969X2

Full Name:

ST6 (alpha-N-acetyl-neuraminy-2,3-beta-galactosyl-1,3)-N-acetylglucosaminide alpha-2,6-sialyltransferase 6

Calculated MW:

333aa, 38 kDa; 374aa, 41 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

