

For Research Use Only

TAF6L Polyclonal antibody

Catalog Number: 15745-1-AP



Basic Information

Catalog Number:

15745-1-AP

Size:

213 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8180

GenBank Accession Number:

BC008785

GeneID (NCBI):

10629

UNIPROT ID:

Q9Y6J9

Full Name:

TAF6-like RNA polymerase II,
p300/CBP-associated factor (PCAF)-
associated factor, 65kDa

Calculated MW:

622 aa, 68 kDa

Observed MW:

68 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB: Jurkat cells, mouse thymus tissue, Raji cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 15745-1-AP (TAF6L antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.