

For Research Use Only

GPN1 Polyclonal antibody

Catalog Number: 15752-1-AP



Basic Information

Catalog Number:

15752-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8252

GenBank Accession Number:

BC007451

GeneID (NCBI):

11321

UNIPROT ID:

Q9HCN4

Full Name:

GPN-loop GTPase 1

Calculated MW:

42 kDa

Observed MW:

42 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HEK-293T cells, human liver tissue, Jurkat cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various cells were subjected to SDS PAGE followed by western blot with 15752-1-AP (GPN1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.